

ERROR FREE I.V. PUMP PROGRAMMING*



SMART I.V. Programmable Pumps

If you are interested in the barcode programmable infusion pumps being offered by Alaris, Sigma Spectrum or Smiths Medical we can assist you. Before you can use a programmable pump you will need a quality label with a 2D symbol that meets the pumps specifications. Dasco works directly with the pump manufacturers to insure your printing system will meet their standards, giving you a 100% guarantee that it will work.


The Dasco solution allows you to connect a printer and generate a quality label that can be scanned to program the pump for drug administration. Reduce errors and increase accuracy with the SMART IV printing system.

CONNECT & PRINT

Intermec printers with the MiSim® software provide a durable, long lasting solution that will interface directly with your pharmacy software. Our integrated smart printers can work with all pharmacy software systems including:

- Misys
- Carecast
- McKesson Star

The Dasco MiSim® software package allows you to connect the printer directly to your current printer port and print labels from your software package. Your label design does not have to change and you can add barcodes. If you would like a quiet printing system that will produce consistent quality labels this is your solution.



Eliminate outdated printers and improve your labeling with a Dasco plug and print system:

- Linear and 2D barcodes
- High Quality printing
- Labels to withstand harsh conditions
- Quiet, high speed printing
- Interface to smart IV pumps

At Dasco, you can choose from our wide variety of specialized stocks to provide you a label that will not smudge or smear and adhere to most any surface. Our systems will allow you to use your current design or modify your label to include barcodes for increased data capture capabilities. If you are having a labeling issue be sure to ask a Dasco specialist for a recommendation to end your headaches.

Dasco has been providing leading edge technology since 1954. Be sure to ask your Dasco representative to develop a solution that meets your needs. For additional information regarding our medical label and printing applications please contact Brad Bellows at 800-927-7701 x105 or via e-mail at bbellows@dascosystems.com.

* The implementation of smart infusion pump systems that utilize drug libraries and 2D imagers can significantly reduce medication administration errors.